

SPECIFICATIONS**CNP-G3 (2-Chloro-4-nitrophenyl maltotrioside)**
AG3-101S

CAS:	118291-90-0
Formula:	$C_{24}H_{34}O_{18}ClN$
Molecular weight:	660.12
Appearance:	Off-white to pale yellow powder
Purity (HPLC):	$\geq 95 \%$
Moisture (K.F):	$\leq 2 \%$
Free CNP:	$\leq 0.1 \%$
Other CNP maltodextrins:	CNPG1 $\leq 0.1 \%$ CNPG2 $\leq 0.1 \%$ CNPG4 $\leq 0.1 \%$
Stability:	At least 1 year at $-20^{\circ}C$