

PRODUCT SPECIFICATIONS

Product :	L-Aspartic acid
Code :	ASP-001
CAS Number :	56-84-8
Formula :	C₄H₇O₄N (Mw: 133.11)
Form :	White powder
Assay :	≥ 99 %
Loss on drying :	≤ 0.2 %
pH (0.5%, H₂O) :	2.5 - 3.5
Heavy metals (Pb) :	≤ 10 ppm
Chlorides (Cl) :	≤ 200 ppm
Specific rotation :	+24.8 - +26.7 °
Stability :	24 months
Storage :	Room Temperature