

PRODUCT SPECIFICATIONS

Product :	Arsenazo III
Code :	ASZ-001
CAS Number :	1668-00-4
Formula :	$C_{22}H_{18}O_{14}As_2N_4S_2$ (Mw: 776.37)
Form :	Deep brown powder
Assay :	> 85 %
Loss on drying :	≤ 8 %
Specific extinction :	300 - 450
Maximum of absorption (in water) :	535 - 543
Maximum of absorption (in buffer, pH 10) :	571 - 576
Calcium :	< 100 ppm
Stability :	24 months
Storage :	Room Temperature