

PRODUCT SPECIFICATIONS

Product :	CAPSO Sodium salt
Code :	CAP-001
CAS Number :	102601-34-3
Formula :	C₉H₁₈O₄NSNa (Mw: 259.3)
Form :	White powder
Assay :	≥ 99 %
Water :	≤ 3 %
pH (1%, diH₂O) :	10.9 – 11.6
A_{280nm} (0.1M H₂O) :	≤ 0.05
A_{260nm} (0.1M H₂O) :	≤ 0.05
Stability :	24 months
Storage :	Room Temperature