

SPECIFICATIONS

Cholesterol esterase from *Candida rugosa* **CHE-210**

Product name	Cholesterol esterase
EC	3.1.1.13.
Appearance	Off white to brown powder
Activity	≥15 U/mg
Optium temperature	37 - 45°C
Optium pH	5.5 - 7.5
Thermal stability	<50°C
Isoelectric point:	4.2
Inhibidor	Sodium azide
Stabilizers	KH ₂ PO ₄
Stability	At least 12 months
Storage	-20 °C

