

PRODUCT SPECIFICATIONS

Product :	Hexokinase
Origin :	Recombinant microorganism
Code :	HK-110
E.C. Number :	2.7.1.1
Nomenclature :	ATP: D-Hexose 6-phosphotransferase
Form :	White to light yellow freeze-dried powder
Activity :	≥ 45 u/mg
Optimum temperature :	45 °C
Optimum pH :	8
Unit definition :	1 unit causes the production of 1 micromole of NADH/min at 37°C and pH 8
Stability :	12 months
Storage :	-20 °C