

PRODUCT SPECIFICATIONS

Product :	D-Lactate Dehydrogenase
Origin :	Microorganism
Code :	LDH-110
E.C. Number :	1.1.1.28
Nomenclature :	(R)-Lactate:NAD⁺ oxidoreductase
Form :	White to greyish powder
Activity :	≥ 200 u/mg
Optimum pH :	6.5 - 8.0
Unit definition :	1 unit causes the oxidation of 1 μmol of NADH/min at 37°C
Stability :	12 months
Storage :	-20 °C