

PRODUCT SPECIFICATIONS

Product :	Magnesium Acetate
Code :	MAT-800
CAS Number :	16674-78-5
Formula :	$\text{Mg}(\text{CH}_3\text{COO})_2 \cdot 4\text{H}_2\text{O}$ (Mw: 214.46)
Form :	White powder
Assay :	$\geq 98 \%$
pH (5%) :	7.5 – 8.5
Sodium :	$\leq 0.5 \%$
Chloride :	$\leq 0.033 \%$
Sulphate :	$\leq 0.06 \%$
Stability :	24 months
Storage :	Room Temperature