

PRODUCT SPECIFICATIONS

Product :	Malate Dehydrogenase
Origin :	Microorganism
Code :	MDH-110
E.C. Number :	1.1.1.37
Nomenclature :	L-Malate:NAD⁺ oxidoreductase
Form :	White to greyish-white powder
Activity :	≥ 40 u/mg
Optimum pH :	8.0 - 8.5
Unit definition :	1 unit causes the oxidation of 1 μmol of NADH/min at 37°C and pH 7.5
Stability :	12 months
Storage :	-20 °C