

PRODUCT SPECIFICATIONS

Product :	Phenol
Code :	PHE-001
CAS Number :	108-95-2
Formula :	C₆H₆O (Mw: 94.11)
Form :	White to light-yellow solid
Assay :	≥ 99 %
Melting point :	38 – 43 °C
Lead (Pb) :	≤ 10 ppm
Cadmium (Cd) :	≤ 1 ppm
Mercury (Hg) :	≤ 1 ppm
Arsenic (As) :	≤ 3 ppm
Stability :	24 months
Storage :	Room Temperature