

PRODUCT SPECIFICATIONS

Product :	p-Nitrophenylphosphate Disodium salt, Hexahydrate (PNPP)
Code :	PPP-001
CAS Number :	4264-83-9
Formula :	$C_6H_4NO_6PNa_2 \cdot 6H_2O$ (Mw: 371.1)
Form :	White to yellow powder
Purity _(HPLC) :	$\geq 98 \%$
Water :	$29 \pm 2 \%$
ϵ_{max}^{312} :	≥ 9600
p-Nitrophenol :	$\leq 0.1 \%$
Free inorganic phosphate :	$\leq 1.0 \%$
Stability :	12 months
Storage :	Room Temperature