

## PRODUCT SPECIFICATIONS

|                             |                                                                |
|-----------------------------|----------------------------------------------------------------|
| <b>Product :</b>            | <b>Acid Cholic, Sodium salt</b>                                |
| <b>Code :</b>               | <b>SCH-001</b>                                                 |
| <b>CAS Number :</b>         | <b>73163-53-8 (alternative: 206986-87-0)</b>                   |
| <b>Formula :</b>            | <b>C<sub>24</sub>H<sub>39</sub>NaO<sub>5</sub> (Mw: 430.6)</b> |
| <b>Form :</b>               | <b>White powder</b>                                            |
| <b>Assay :</b>              | <b>≥ 98 %</b>                                                  |
| <b>pH (2%) :</b>            | <b>7.5 – 8.5</b>                                               |
| <b>Sodium deoxycholate:</b> | <b>≤ 2 %</b>                                                   |
| <b>Loss on drying :</b>     | <b>≤ 5 %</b>                                                   |
| <b>Stability :</b>          | <b>24 months</b>                                               |
| <b>Storage :</b>            | <b>Room Temperature</b>                                        |