

PRODUCT SPECIFICATIONS

Product :	Sodium Taurodeoxycholate
Code :	STD-001
CAS Number :	207737-97-1
Formula :	C₂₆H₄₄NO₆SNa (Mw: 521.7)
Form :	White powder
Assay :	≥ 97 %
Loss on drying :	≤ 7.0 %
Sodium deoxycholate :	≤ 1.0 %
Sodium taurocholate :	≤ 0.5 %
Taurine :	≤ 1.0 %
Bile acid sodium :	≤ 0.5 %
Heavy metals :	≤ 20 ppm
Stability :	24 months
Storage :	Room Temperature